

NU-DERM™
ALGINATE WOUND DRESSING

Designed for
high absorption



WHAT IS IT?

NU-DERM™ Alginate Wound Dressing is a sterile, non-woven pad consisting of high G (guluronic acid) alginate and carboxymethylcellulose (CMC) fiber. NU-DERM™ Alginate Wound Dressing is designed to be fast gelling and to allow intact removal.

HOW DOES IT WORK?

High absorbency

NU-DERM™ Alginate is a highly absorbent dressing. *In vitro* data shows that it can absorb up to 18 times its own weight in fluid, which may reduce the frequency of dressing changes.


NU-DERM™ Alginate forms a hydrophilic gel at the wound surface to help create a moist wound environment.



The average* *in vitro* absorbency¹ for NU-DERM™ Alginate is 23.9 g/100cm²

Wet tensile strength

In vitro data demonstrates that dressing tensile strength is maintained when wet¹ to facilitate intact removal and minimize the need for wound irrigation at dressing change.



Intact dressing removal

WHEN TO USE


NU-DERM™ Alginate Wound Dressing is indicated for the management of moderate to heavily exuding chronic wounds and to control minor bleeding in superficial acute wounds.

Chronic wounds: pressure ulcers, leg ulcers, venous stasis ulcers, diabetic ulcers and arterial ulcers.

Superficial acute wounds: abrasions, lacerations, donor sites and postoperative wounds.

NU-DERM™ Alginate Wound Dressing can be cut or folded to accommodate the wound bed.

Item code	Size	Systagenix unit of measure	HCPCS code
AWD112	1" x 12"	5 ea/ct, 5 ct/bx	A6199
AWD202	2" x 2"	10 ea/ct, 5 ct/bx	A6196
AWD404	4" x 4"	10 ea/ct, 5 ct/bx	A6196
AWD408	4" x 8"	5 ea/ct, 5 ct/bx	A6197



To learn more about the benefits of NU-DERM™ Alginate contact your Acelity representative at **800-275-4524** or visit acelity.com