

NU-DERM™
HYDROCOLLOID WOUND DRESSING

A hydrocolloid in a range
of **shapes** and **sizes** for
your patients.



WHAT ARE THEY?

NU-DERM™ Hydrocolloid Dressings are designed to help maintain a moist wound environment. There are three variants in the range, available in a variety of shapes and sizes. All dressings are waterproof and will remain in place during showering.

NU-DERM™ Hydrocolloid Dressing Range:

NU-DERM™ Border Hydrocolloid Wound Dressing

It is conformable, has a top layer of low-friction film, and has bevelled edges all around the product - the border itself is a continuation of the skin-friendly hydrocolloid adhesive material.

NU-DERM™ Standard Hydrocolloid Wound Dressing

It has a top layer of foam, rounded corners, and is uniform in thickness.

NU-DERM™ Thin Hydrocolloid Wound Dressing

The semi-transparent variance, which is conformable and has a top layer of low-friction film, allows close and easy monitoring of the wound bed without dressing removal.

HOW DO THEY WORK?

NU-DERM™ Hydrocolloid Dressings help maintain a physiologically moist wound environment conducive to wound healing.

NU-DERM™ Hydrocolloid Dressings helps protect wounds against external contamination.

2x

Fluid handling

In vitro data show that NU-DERM™ Hydrocolloid Dressings fluid handling capacity in a period of 48h is twice as much as DuoDERM® CGF® Sterile Dressing and DuoDERM Signal® Dressing in (NU-DERM™ Hydrocolloid Dressings 6.24g/48h, DuoDERM® Dressing 2.31, DuoDERM Signal® Dressing 2.91g/48h)¹.

WHEN TO USE THEM?

NU-DERM™ Border Hydrocolloid Dressing and NU-DERM™ Standard Hydrocolloid Dressing are primarily indicated for the management of lightly to moderately exuding pressure sores and leg ulcers.

NU-DERM™ Thin Hydrocolloid Dressing is primarily indicated for the management of superficial dry to lightly exuding wounds, post-operative wounds, and superficial wounds and abrasions; however it is also useful on small wounds towards the end of the healing phase.

NU-DERM™ Hydrocolloid Dressings			
Item code	Size	Eaches/UOM	HCPCS code
NU-DERM™ Thin Hydrocolloid Dressing			
HCT101	4" x 4"	10 ea/tray, 10 trays/ca	A6234
NU-DERM™ Standard Hydrocolloid Dressing			
HCF204	4" x 4"	5 ea/tray, 10 trays/ca	A6234
HCF204	8" x 8"	5 ea/tray, 4 trays/ca	A6236
NU-DERM™ Border Hydrocolloid Dressing			
HCB102	2" x 2"	20 ea/tray, 5 trays/ca	A6237
HCB204	4" x 4"	10 ea/tray, 16 trays/ca	A6237
HCB204	6" x 6"	5 ea/tray, 4 trays/ca	A6238
NU-DERM™ Sacral Hydrocolloid Dressing			
HCS100	6" x 7"	5 ea/tray, 4 trays/ca	A6238

To learn more about the benefits of NU-DERM™ Hydrocolloid Dressings contact your Acelity representative at **800-275-4524** or visit **acelity.com**

NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for Systagenix products. Please consult a healthcare provider and product instructions for use prior to application. Rx only.

Reference

1. Data on file

Copyright 2016, 2017 Systagenix Wound Management, Limited. All rights reserved. DuoDerm is a trademark of Convatec. All other trademarks designated herein are proprietary to Systagenix Wound Management, Limited, its affiliates and/or licensors. LIT#NUDERM HYD SS • PRA000593-RO-US, EN (05/17)



Acelity™