

Onli

Ready to Use Hydrophilic Intermittent Catheter

Whether their day is ordinary or extraordinary, the Onli catheter is there for them

Designed to easily fit into your patient's daily routine as a hassle-free and mess-free experience

Ready. Simple.





Onli

Ready to Use Hydrophilic Intermittent Catheter

Relief they can count on, whatever the day may bring

The Onli catheter is designed for easy, everyday use


- The Onli ready to use hydrophilic intermittent catheter is designed to provide a hassle-free and mess-free experience
- As easy to learn as it is to use, you can be confident that your patients will be able to use the Onli catheter the right way at the right time



The catheter design provides a balance of flexibility and rigidity to support ease of insertion

- 1 Prehydrated and ready to use** with reduced spill¹ and no stain^{*2}
- 2 Smooth hydrophilic coating** for easy insertion and withdrawal
- 3 Ergonomic gripper** on 16" catheter provides better control
- 4 Easy open package**

^{*}on denim, linen and polyester



Onli

Ready to Use Hydrophilic Intermittent Catheter

Onli Ready to Use Hydrophilic Intermittent Catheter

Stock No	Length (inches)	Size	Box Qty
82081-30	7	8Fr	30
82101-30	7	10Fr	30
82121-30	7	12Fr	30
82141-30	7	14Fr	30
82084-30	16	8Fr	30
82104-30	16	10Fr	30
82124-30	16	12Fr	30
82144-30	16	14Fr	30
82164-30	16	16Fr	30



Call today for more information or a sample: 888.740.8999 or visit hollister.com/onli

This intermittent catheter is a flexible tubular device that is inserted through the urethra by male, female and pediatric patients who need to drain urine from the bladder.

- Ready to use, single-use intermittent catheter
- Hydrophilic-coated, phthalate-free catheter
- 2 smooth catheter eyelets
- Color-coded funnel
- Not made with natural rubber latex
- Gripper on 16-inch catheter to aid insertion

WARNING: To help reduce the potential for infection and/or other complications, do not reuse.

WARNING: If discomfort or any sign of trauma occurs, discontinue use immediately and consult your healthcare practitioner.

PRECAUTION: Please consult your healthcare practitioner before using this product if any of the following conditions are present: • Severed urethra • Unexplained urethral bleeding • Pronounced stricture • False passage • Urethritis – inflammation of the urethra • Prostatitis – inflammation of the prostate gland • Epididymitis – inflammation of the epididymis (testicle tube)

PRECAUTION: Self-catheterization should follow the plan of care and advice given by your healthcare practitioner and be carried out only in accordance with the instructions for use provided. Because catheterization frequency varies by person, the recommended frequency of your catheterization should be provided by your healthcare practitioner. For any other questions about your catheterization, please contact your healthcare practitioner.

CAUTION: Rx Only. Federal (USA) Law restricts this device to sale by or on the order of a physician or other healthcare practitioner licensed under state law to order this product. Prior to use of this device, be sure to read (i) the complete information on how to use this device including Warnings, Cautions, and Instructions for Use, and (ii) all other packaging inserts and labels supplied with the product.

NOTE: Store boxes in a flat position and at normal room temperature.

This product consists of: thermoplastic gripper (16-inch catheter only), polyvinylchloride funnel, and tubing not made with phthalates.

1. Data on file. Loose water in Onli catheters and competitor products packaging study. Hollister, 2015.
2. Data on file. Analysis of potential loose liquid to stain clothing for Onli catheters and competitor products. Hollister, 2015.



Hollister Incorporated
2000 Hollister Drive
Libertyville, Illinois 60048
1.800.323.4060

www.hollister.com

The Hollister logo and Onli are trademarks of Hollister Incorporated
© 2017 Hollister Incorporated. 923093-117